

Wherever you're going, At'm goes, too. It's the power chair lightweight and compact enough to fit in the trunk of a car. Simple enough to assemble in 60 seconds. Powerful enough for a full day of adventure.

How easy is At'm? If you can set up a lawn chair, you're ready. No mobility lifts, no tools, just three easy pieces to assemble. Open the seat, snap it on the lightweight base, add the battery and go! No wonder we call At'm the Take Along Chair!



Invacare® At'm Take Along Chair

The At'm goes wherever you go. Folds together like a lawn chair. Goes together in seconds. How easy is it to drive? Flip the switch and go.



SPECIFICATIONS

Dimensions	
Overall width	25.5"
Overall length	35.5"
Seat style	Sling style
Seat depth	14", 16"
Seat width	16", 18"
Seat back height	18", 21"

Tires	
Front tire size	1" x 6"
Rear tire size	2.25" x 12.5"

Weight	
Base section	34 lb.
Seat and arms	14 lb.
Battery pack	26 lb.
Total weight	74 lb.
Weight limit	250 lb.

Mechanical/Performance	
Speed range	0 - 4 mph
Incline capability	6°
Ground clearance	2"
Charger	off-board
Battery size	17 Ah
Battery range	10 - 12 miles*

Limited Warranty	
Frame	36 months
Electronics	12 months

*The Invacare® At'm is not intended for street or highway use. MINI is a registered trademark of Bayerische Motoren Werke Aktiengesellschaft Corporation and is not affiliated with Invacare Corporation.

Available colors



Candy Red Pearl



Deep Blue Pearl



Take Along Chair



Yes, you can.™

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare and the blue medallion design are registered trademarks and Yes, you can. is a trademark of Invacare Corporation. Specifications are subject to change without notification.

Patent pending.

© 2004

Invacare Corporation
Form No. 04-159



"Invacare says 'Yes, you can™,' and I couldn't have said it better myself."

